


## IKA Vibrax VXR Basic



wide speed range (0 - 2200 rpm)

Rating: Not Rated Yet

**Price**

Sales price \$550.00

[Ask a question about this product](#)

Manufacturer [IKA](#)

### Description

Optoelectronically controlled small shaker with a very wide speed range (0 - 2200 rpm).

### Features:

- Suitable for continuous operation
- New design and improved drive system
- Circular shaking motions
- Slow speeds are well maintained

## Shakers/Stirrers/Hotplates: IKA Vibrax VXR Basic

---

Attachments are interchangeable

Specifications:

Type of movement: orbital  
Shaking stroke: 4 mm  
Permissible shaking weight (incl. attachment): 2 kg  
Motor rating input: 35 W  
Motor rating output: 13.5 W  
Permissible ON time: 100 %  
Speed min (adjustable): 100 rpm  
Speed range: 0 - 2200 rpm  
Speed display: scale  
Operating mode: continuous operation  
Dimensions (W x H x D): 157 x 130 x 247 mm  
Weight: 6.1 kg  
Permissible ambient temperature: 5 - 40 °C  
Permissible relative humidity: 80 %  
Voltage: 100 - 120 V  
Frequency: 50/60 Hz  
Power input: 45 W

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)