

## IKA Works Ceramag Midi



Speed: 150-1200 RPM

Rating: Not Rated Yet

**Price**

Sales price \$225.00

[Ask a question about this product](#)

Manufacturer [IKA](#)

### Description

Hotplate stirrer that provides you with a flexible and carefully monitored heating and stirring. Independent controls for heat and stirring, for individual needs or simultaneous use.

### Specifications:

Speed: 150-1200 RPM

Plate Dimensions: 5.5 x 5.5

Temperature Range: 0 to 450C

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)

