

Infors MultitronTriplestack



Incubator shaker w/humidity, CO2, and 25mm orbit.

Rating: Not Rated Yet

Price

Sales price \$27,500.00

[Ask a question about this product](#)

Description

Flexibility combined with user-friendly operation and operational safety with the optimum utilization of space. Easy-to-use: tray handling, auto-restart, front door opening, and clear operating panel. Can hold up to 5L flasks. Even temperature distribution over the whole shaking tray. Safe, stable stacking of three incubator shakers at one time.

Incubator Shakers: Infors MultitronTriplestack

Specifications:

Temperature Range: 5 above ambient to 60C
Shaker Speed: 30x400 RPM
Platform Size: 840 x 450mm
Maximum Capacity: 6 x 5000ml Flasks
Orbital Diameter: 25 - 50 mm
Dimensions (W x D x H): 1070 x 860 x 550 mm

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com