

Integra Mini 96 Pipettor



Designed for 96- and 384-well microplates

Rating: Not Rated Yet

Price

Sales price \$7,750.00

[Ask a question about this product](#)

Description

The INTEGRA Mini 96 Pipettor is a lightweight, handheld 96-channel electronic pipette that enables rapid, accurate filling of 96- and 384-well plates, boosting productivity for high-throughput workflows.

Specifications

Channels: 96 (air-displacement)

Volume Ranges: Four models: 0.5–12.5 μL , 5–125 μL , 10–300 μL , or 50–1250 μL

Accuracy & Precision: Accuracy 1.2–1.8 %, precision \approx 1.0 % (varies by model)

Portability: \approx 9.4 kg with carry handles; compact footprint fits laminar flow cabinets

Plate Compatibility: Designed for 96- and 384-well microplates; supports full/partial plate filling with optional two-position stage

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office

650-276-7487 fax
sales@biotechequipmentsales.com