

## Invitrogen PowerEase Touch 600



4.3-inch backlit LCD with pre-programmed protocols

Rating: Not Rated Yet

**Price**

Sales price \$1,150.00

[Ask a question about this product](#)

### Description

The Invitrogen PowerEase Touch 600W (models PS0600 / PS0601 etc.) is a high-power electrophoresis power supply from Thermo Fisher Scientific. It features a modern 4.3-inch color touchscreen interface and delivers up to 500 V, 3 A (3000 mA), and 600 W. This makes it ideal for demanding applications such as running and transferring mini/midi protein gels (SDS-PAGE, native PAGE), IEF, western blotting, and DNA/RNA electrophoresis. It supports simultaneous operation of multiple gel tanks via four output pairs. This is a powerful, user-friendly unit well-suited for busy molecular biology labs.

### Key Features & Specifications

**High Output Capacity** — Maximum 500 V, 3 A (3000 mA), and 600 W; supports constant voltage, constant current, or constant power modes with programmable limits

**Intuitive Touchscreen Interface** — 4.3-inch backlit LCD with pre-programmed protocols for common gels/transfers, custom method creation, timers (up to 999 min), and real-time monitoring

**Four Output Pairs** — Four sets of terminals (4 positive + 4 negative) allow running/transferring multiple electrophoresis units simultaneously.

**Safety & Usability Features** — Over-voltage, over-current, over-temperature, and no-load protection; stackable design with adjustable viewing angle feet; audible alarms and user-adjustable settings (brightness, beep volume, etc.).

**Compact & Versatile Design** — Dimensions approx. 9 x 12 x 4 inches (23 x 30 x 10 cm); weight ~2.8 kg; universal power input (115/230 VAC, 50/60 Hz); suitable for both routine and high-current applications like midi gels and blotting.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc. has been based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.

South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)