


## Beckman J-LITE®JLA-16.250 Fixed-Angle Aluminum Rotor, Single Lock-6x250mL, 16,000 rpm, 38,400 x g



Lightweight alternative to conventional 250 rotors

Rating: Not Rated Yet

### Price

Sales price \$2,650.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

### Description

JLA-16.250 Fixed-Angle Rotor with Single Lock Lid, Aluminum.

General purpose, large-volume and multitube processing. Lightweight alternative to conventional 250 rotors; empty rotor weight is 10.3 kg (22.7 lb). For use in Avanti® J Series and J2 Series centrifuges. The maximum speed is 14,000 rpm in J2 Series centrifuges.

#### Specifications:

Platform: Avanti

Rotor Type: Fixed-Angle (High Performance)

Maximum RPM: 16,000 rpm

k-Factor: 1,090

Nominal Capacity: 1,500.00|mL

rMAX: 134

rMIN: 46

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)