

## KD Scientific Model 210



KD Scientific Modle 210 infusion/withdrawal dual syringe pump is BRAND NEW.

Rating: Not Rated Yet

**Price**

Sales price \$1,650.00

[Ask a question about this product](#)

Manufacturer [KD Scientific](#)

### Description

**!!!SOLD!!!**

This KD Scientific Modle 210 infusion/withdrawal syringe pump is BRAND NEW. It's microprocessor controlled, hold two syringes, has a broad speed range, holds a wide range of sizes, has permanent memory of last setting, continuous dispense volume display, includes an optical encoder stall detector, built-in RS-232C interface, TTL interface, has 5 operation modes, with independent rate & volume settings for both infusion & withdrawal.

0.001 microliters/hour minimum flow rate w/10 microliter syringe, 21 ml/hr maximum flow rate w/10 ml syringe - 71 ml/min w/60 ml syringe

1600 steps/sec maximum rate

**Please call us at (650) 871-5707 to check availability of this product!**