

KNF Lab Laboport Vacuum Pump



Delivery rate at atm. pressure [l/min]* max.: 34

Rating: Not Rated Yet

Price

Sales price \$2,250.00

[Ask a question about this product](#)

Description

Model UN842.3FTP

Specifications:

Pneumatic performance

Max. permissible operating pressure [bar g]: 1.0

Ultimate vacuum [mbar abs.]: 2

Delivery rate at atm. pressure [l/min]* max.: 34

Pneumatic connection: Hose connection [mm] ID 10

Ambient and media temperature:

Permissible ambient temp.: + 5 °C to + 40 °C

Permissible media temperature: + 5 °C to + 40 °C

Other parameters:

Weight [kg]: 13.4

Dimensions: L x H x W [mm] 341 x 223 x 167

Maximum permissible ambient relative humidity: 80 % for temperatures up to 31 °C, decreasing linearly to 50 % at 40 °C

Maximum altitude of site [m above sea level]: 2000

115V / 60 Hz / 3.2 A / 290 W

Fuse T6.3A

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com