

Lab-Line Ultra-Clean 100



Lab-Line Ultra-Clean 100 Cleanroom Oven, mechanical convection, microprocessor-controlled, 4.6 cu/ft

Rating: Not Rated Yet

Price

Sales price \$2,750.00

[Ask a question about this product](#)

Manufacturer [Lab-Line](#)

Description

!!!SOLD!!!

This model 3499M-3 Lab-Line Ultra-Clean 100 Cleanroom Oven maintains Class 100 conditions without the use of filters at operating temperatures up to 250C, all stainless interior, 2 built-in safety systems to protect samples, visual and audible overtemp alarms microprocessor-controlled, 4.6 cu/ft, 208V, 50/60 Hz, 3100W