

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

LabArmor 20L Bead Bath w/out Cover, Model 74300-720



Capacity: 20 L (0.71 cu. ft.); +5°C to 80°C

Rating: Not Rated Yet

Price

Sales price \$700.00

[Ask a question about this product](#)

Description

LabArmor 20L Bead Bath (model 74300-720, without cover) is a larger-capacity waterless digital laboratory bath that uses thermal beads (sold separately) in place of water. It delivers superior temperature uniformity, eliminates contamination risks from water baths, reduces energy consumption (up to 50% less), and requires virtually no maintenance—no refilling, evaporation, or germicides.

Specifications

Capacity: 20 L (0.71 cu. ft.); recommended bead volume: ~16 L

Temperature Range: Ambient +5°C to 80°C; resolution 0.1°C; stability $\pm 0.1^{\circ}\text{C}$; uniformity $\pm 1.0^{\circ}\text{C}$

Digital Control: Microprocessor-based with clear display for precise, stable operation

Dimensions: Exterior: 15.9" W x 23.4" D x 8.5" H (404 x 594 x 216 mm); Interior: 12" W x 17.1" D x 6" H.

Electrical & Weight: 120V / ~6.0A (50/60 Hz); unit weight 22.0 lbs (10.0 kg) without beads. Stainless steel interior, painted steel exterior

This model suits higher-volume warming, thawing, or incubation needs in labs. Note: Beads are required for proper function and not included; a

cover is available separately.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com