

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

LabArmor 6L Bead Bath w/ Cover, Model 74300-706



Capacity: 6 L (0.21 cu. ft.); Temperature Range: Ambient +5°C to 80°C; resolution 0.1°C; stability $\pm 0.1^\circ\text{C}$ (at 37°C); uniformity $\pm 1.0^\circ\text{C}$

Rating: Not Rated Yet

Price

Sales price \$500.00

[Ask a question about this product](#)

Manufacturer [Shel Lab](#)

Description

LabArmor 6L Bead Bath (model 74300-706) is a compact, waterless digital laboratory bath that uses thermal beads (sold separately) instead of water. It provides excellent temperature uniformity, eliminates microbial contamination risks associated with traditional water baths, reduces energy use (up to 50% less), and requires minimal maintenance—no evaporation, refilling, or germicides.

Specifications

Capacity: 6 L (0.21 cu. ft.); recommended bead volume: 5 L

Temperature Range: Ambient +5°C to 80°C; resolution 0.1°C; stability $\pm 0.1^\circ\text{C}$ (at 37°C); uniformity $\pm 1.0^\circ\text{C}$

Digital Control: Microprocessor-based with easy-to-read display; includes stainless steel lid/cover

Dimensions: Exterior: 15.9" W x 12.3" D x 8.5" H (404 x 312 x 216 mm); Interior: 12" W x 6.1" D x 6" H

Electrical: 120V / 3.0A (50/60 Hz); unit weight ~13 lbs (5.9 kg) without beads.

Baths/Circulators: LabArmor 6L Bead Bath w/ Cover, Model 74300-706

This model is ideal for warming, incubating, or thawing samples in a clean, efficient manner. Note: Requires LabArmor beads for proper operation (not included).

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com