

Labconco Purifier Delta Series



Labconco 5 foot Purifier Delta Series biosafety cabinet on stand

Rating: Not Rated Yet

Price

Sales price \$3,500.00

[Ask a question about this product](#)

Manufacturer [Labconco](#)

Description

!!!SOLD!!!

SOLD!!!!

This Labconco 5 foot Purifier Delta Series biosafety cabinet comes w/its stand and is in very good condition. It's a Class II, Type A2 & suitable for Biosafety Levels 1, 2, & 3 as well as ISO Class 5. All stainless interior.

Please call us at (650) 871-5707 to check availability of this product!