

## LAUDA-Brinkmann WKL 230 Circulation Chiller



Circulation chiller with closed cooler circuit for continuous nonstop work in research, technology and production at temperatures from -10 up to 40 °C.

Rating: Not Rated Yet

**Price**

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Lauda-Brinkmann](#)

### Description

LAUDA WK class circulation chillers are used where operating heat has to be discharged reliably and quickly in chemical production processes or from technical systems – in harsh constant use environments. Unlike cooling with tap water, they provide a constant temperature irrespective of the time of year and pressure fluctuations. The water saving also protects the environment and significantly reduces operating costs. The WKL 230 circulation chiller was especially designed for simple thermostating tasks below 0 °C. Thanks to its extremely compact construction and small footprint it fits anywhere. The temperature range of the WKL 230 is from -10 °C up to 40 °C with a temperature stability of  $\pm 0.5$  K.