

Leica SM 2400 Microtome



Microtome with safety features.

Rating: Not Rated Yet

Price

Sales price \$4,250.00

[Ask a question about this product](#)

Manufacturer [Leica](#)

Description

Leica SM2400 is an instrument that can be used for a variety of applications in all laboratories where paraffin embedded materials need to be sectioned.

The safety features of the system include a sledge locking device and knife guard.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com