

## Lindberg/Blue M Mini-Mite



Lindberg/Blue M Mini-Mite tube furnace to 1100C, 12" x 1" heated length, lightly used

Rating: Not Rated Yet

**Price**

Sales price \$1,050.00

[Ask a question about this product](#)

Manufacturer [Lindberg](#)

### Description

!!!SOLD!!!

This Thermo Lindberg/Blue M TF55030A Mini-Mite tube furnace gets to 1100C has a 12" x 1" heated length, is controlled by a single setpoint controller, has an easy to use split hinge design, and is lightly used w/a small dent on the top rear of the instrument (that doesn't affect function at all).