

Lindberg Blue M Oven



Temperature range: 40-300C

Rating: Not Rated Yet

Price

Sales price \$1,400.00

[Ask a question about this product](#)

Manufacturer [Lindberg](#)

Description

Model: MO1440A-1

Features internal blowers, 2 stainless steel shelves.

Specifications:

Temperature range: 40-300C

Temperature uniformity: 1% at +100C

Heat-up Rate: 15 min to 100C; 44 min to 200C; 115 min to 300C

Interior Dimensions (W x D x H): 22 x 16 x 15 in.

Exterior Dimensions (W x D x H): 34 x 22.5 x 37 in.

Electrical requirements: 120V / 50/60Hz / 1900W

Ovens/Furnaces: Lindberg Blue M Oven

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com