

Mettler Toldeo PG5002-SDR



Mettler Toldeo PG5002-S Standard Level toploading balance, Delta Range, 5100g/1000g x 0.1/0.01g

Rating: Not Rated Yet

Price

Sales price \$1,050.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

!!!SOLD!!!

This Mettler Toldeo PG5002-S Standard Level toploading balance, Delta Range, 5100g/1000g x 0.1/0.01g offers the same weighing technology and rugged construction as the Professional Level balances. It features FACT (Fully Automatic Calibration Technology) which delivers automatic calibration whenever ambient temperature changes by +/- 2C, has a built in RS232 interface, dynamic weighing, full range taring, large 6.5" x 6.5" weighing pan, and multiple weighing modes.

Please call us at (650) 871-5707 to check availability of this product!