

Mettler Toledo AG204DR



Mettler Toldedo AG204DR Delta Range Analytical balance, 61g/210g by 0.1mg/1mg

Rating: Not Rated Yet

Price

Sales price \$1,325.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

!!!SOLD!!!

This Mettler Toldedo AG204DR Delta Range Analytical balance, 61g/210g by 0.1mg/1mg has Fully Automatic Calibration Technology (FACT), built-in universal interface, a 200g internal calibration weight, and an 8.5 cm diameter weighing pan.

Comes with full range taring. Units of measure include: g, mg, oz, oz t, GN, dwt, ct, S tl, H t, & T tl

Please call us at (650) 871-5707 to check availability of this product!