

## Mettler Toledo Balance Printer LC-P45



LC-P45 printer for Mettler Toledo balances and other products.

Rating: Not Rated Yet

**Price**

Sales price \$500.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

### Description

The LC-P45 is a versatile dot matrix printer with advanced print functions for normal paper. It can supplement the weight value with date, time and other data. It has two LocalCAN interfaces and one RS-232-C interface and is primarily suitable for attachment to Mettler Toledo balances and scales equipped with the standard software, but will also print out the records of the Mettler Toledo application software Pacs (M-Pacs) and of the LP16 and LJ16 Mettler Toledo Infrared Dryers. The various print functions of the LC-P45 can be selected in the dialog mode (English, German, French, Spanish or Italian). The settings are stored in a failsafe manner thanks to battery buffering.

[Biotech Equipment Sales, Inc.](#)  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)