

Mettler Toledo ME104TE Analytical Balance



120g x 0.1 mg

Rating: Not Rated Yet

Price

Sales price \$1,550.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

Rugged, touchscreen-controlled analytical balance offering high-precision weighing up to 120 g with 0.1 mg readability.

Key Specifications

Capacity & Readability: 120 g maximum capacity with 0.1 mg (0.0001 g) readability

Repeatability: Typical ± 0.08 mg; guaranteed ± 0.1 mg

Linearity: ± 0.2 mg maximum, typically ± 0.06 mg

Minimum Weight (USP, 0.1 %, k=2): 160 mg

Pan Size: 90 mm diameter stainless-steel pan

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-

quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com