

Mettler Toledo ME2002E Top-loading Balance



Capacity x Readability: 2200g x 0.01g

Rating: Not Rated Yet

Price

Sales price \$895.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

Top-loading balance with a bright display providing easy access to settings, features, and applications, providing everything you need for fast and efficient weighing. The solid metal base, stainless-steel weighing pan and overload protection ensure reliable performance, even under intense daily use.

2200 g capacity; 10 mg readability; brilliant display; 1 x RS232; external adjustment; solid metal base; easy levelling; 10 applications.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080

650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com