

Mettler Toledo ML1602T Precision Balance



Maximum capacity: 1620 g

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

The Mettler Toledo ML1602T (often called 1602T) is a compact ML-T precision balance for routine lab weighing. It uses a MonoBloc cell with FACT automatic internal adjustment, a 4.5-inch color touchscreen, user management, and built-in applications such as counting, check weighing, and percent weighing.

Key specifications

Maximum capacity: 1620 g

Readability: 0.01 g

Settling time: 1.5 s (typical repeatability 0.007 g; linearity ± 0.02 g)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com