


Mettler Toledo MS104TS/00 Analytical Balance



Maximum capacity: 120 g

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

The Mettler Toledo MS104TS/00 is an advanced-level analytical balance in the MS-TS line (material no. 30133522). It uses a MonoBloc weighing cell with automatic internal adjustment (FACT), a 7-inch color touchscreen, user management, and 11 built-in applications. It is aimed at chemical, pharma, and food & beverage labs.

Key specifications

Maximum capacity: 120 g

Readability: 0.1 mg

Settling time: 2 s (typical repeatability 0.08 mg; linearity ± 0.2 mg)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com

