

## Mettler Toledo Portable Lab Densito 30PX New in Case



Quick, simple measurements; 30PX / 30GS New in Case

Rating: Not Rated Yet

**Price**

Sales price \$950.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

### Description

Portable, handheld measuring instrument analyzes liquids for density or specific gravity. Mettler Toledo Densito 30PX Density Meter operates on the oscillating tube principle. Results are automatically calculated into one of the following units: density, specific gravity, API degrees, Brix, % Alcohol, % H2SO4, "Baume," Plato, Proof or a user defined unit. Shows results within seconds.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)