


Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Mettler Toledo Precision Balance ME4001E



4200 g maximum capacity with 0.1 g (100 mg) readability

Rating: Not Rated Yet

Price

Sales price \$950.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

Mettler Toledo ME4001E is a robust, entry-level precision balance from the ME series, designed for reliable everyday weighing in laboratories, education, and industrial settings. It offers essential functionality with a solid metal base for durability and straightforward operation.

Key Features:

Capacity & Readability: 4200 g maximum capacity with 0.1 g (100 mg) readability.

Display & Operation: Bright, clear LCD display with simple menu navigation and basic operation keys for fast, intuitive use.

Adjustment: External calibration/adjustment for accuracy and traceability.

Connectivity: 1 x RS232 interface for data transfer to printers or computers.

Design: Solid metal base, stainless steel weighing pan, overload protection, and level indicator for stable performance. Supports applications

Balances: Mettler Toledo Precision Balance ME4001E

like weighing, piece counting, percent weighing, and more.

This balance emphasizes reliability, ease of use, and value for routine tasks.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com