

## Mettler Toledo SB32001 DeltaRange Top-loading Precision Balance



0.1 g / 1 g readability and high capacity (6400g / 32,100g)

Rating: Not Rated Yet

**Price**

Sales price \$1,950.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

### Description

Designed for weighing in laboratory or industrial environment with 0.1 g / 1 g readability and high capacity (6400g / 32,100g)

Full metal housing and a high chemical resistance. Simple to operate. MonoBloc weighing cell.

The balance features moveable fine range with ten times better readability as well as large, stainless steel weighing platform.

Size of weighing pan:349x232

Dimensions:381x321x92 mm (WxDxH)

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)