

Mettler Toledo SR16000



Mettler Toldedo SR16000 Professional Level top-loading balance, 16,100g x 1g

Rating: Not Rated Yet

Price

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

!!!! - RECENTLY SOLD - !!!!

This Mettler Toldedo SR16000 Professional Level top-loading balance, 16,100g x 1g, has an 11" x 14.25" weighing pan, a SmartBar, which provides quick access to all balance functions, ProFACT internal calibration w/history function, initiated by 1.5C temp. fluctuation, and full range taring