

Mettler Toledo XP2003SDR Mass Comparator



Capacity: 2,300 g; Readability: 1 mg

Rating: Not Rated Yet

Price

Sales price \$1,250.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

High-precision mass comparator.

Technical Specifications:

Capacity: 2,300 g

Readability: 1 mg

Pan Size: 125 x 125 mm

Repeatability (standard deviation of 10 comparison weighing ABA after elimination of drift): 1 mg

Typical repeatability: 0.8 mg + $9 \times 10^{-8} \sqrt{C \times R_{gr}}$

Balances: Mettler Toledo XP2003SDR Mass Comparator

Linearity: 5 mg
Typical stabilization time: 2?Ç³10 s
Adjustment: proFACT
Adjustment with external weight: >=500 g
Sensitivity drift±10-6
Free space above weighing pan: 235 mm
Feedthrough for below-the-balance weighing: Standard
Backlit color graphics display with touch screen operation: Standard
WeighCom application for guided mass comparison application: Standard
Temperature (°C): 10?Ç³30 °C
Max. temperature change: (°/24 h)0.5
Relative air humidity (%): (±2%)40-70
Dimensions (WxDxH) / Net weight: 214 x 257 x 363 mm / 7.1. kg
Display unit: 194 x 133 x 58 mm / 1.2 kg
AC adapter: 108 x 67 x 37 mm / 0.3 kg

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com