

Mettler Toledo XPR404S Analytical Balance



Capacity of 410 g; readability of 0.0001 g; Ø2.7 (68 mm) pan size.

Rating: Not Rated Yet

Price

Sales price \$5,350.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

Mettler-Toledo XPR small platform analytical balances deliver outstanding weighing performance and support the highest requirements for data integrity and compliance.

Technical Specifications:

Maximum Capacity: 410 g

Readability: 0.1 mg

Repeatability (typical): 0.06 mg

Minimum Weight (USP), Typical: 120 mg

Weighing Pan Dimensions (DxW): 68 x 68 mm (SmartPan)

Balances: Mettler Toledo XPR404S Analytical Balance

Linearity: ± 0.2 mg
Repeatability (Test Weight): 0.1 mg (400 g)
Settling Time: 2 seconds
Dimensions (DxHxW): 411 x 368 x 214 mm
Adjustment: Internal / proFACT Advanced
Linearity: 0.2 mg
Material Number(s): 30315628

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.

**226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com**