

Balances: Mettler Toledo XPR6002S w/P-52RUE Printer

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Mettler Toledo XPR6002S w/P-52RUE Printer



6100 g (6.1 kg) capacity with 0.01 g (10 mg) readability

Rating: Not Rated Yet

Price

Sales price \$4,650.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

Mettler Toledo XPR6002S with P-52RUE Printer is a high-performance precision balance from the XPR Excellence series, paired with a reliable dot matrix data printer. This combination delivers accurate, traceable weighing results with seamless documentation for regulated labs, quality control, and research applications.

Key Features:

Capacity & Readability: 6100 g (6.1 kg) capacity with 0.01 g (10 mg) readability; supports large samples with excellent accuracy.

SmartPan & Performance: SmartPan weighing pan for faster, more stable results and error reduction; internal calibration for consistent performance.

Intuitive Interface: 7-inch color capacitive touchscreen with glove-compatible operation, built-in methods, results notepad, and user management/passwords.

Connectivity: Multiple ports (4x USB, 1x LAN/Ethernet) for error-free data transfer, plus full integration with the P-52RUE printer.

P-52RUE Printer: Compact dot matrix printer for direct output of weighing results, labels, reports, and GLP/GMP-compliant records via

RS-232/USB/Ethernet; reliable long-term printing at up to 10 mm/s.

This setup excels in compliance, efficiency, and ease of use for routine or complex weighing tasks.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com