

Mettler Toledo XSR104 Excellence Analytical Balance



Maximum capacity: 120 g

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

The Mettler Toledo XSR104 is a lab balance with a hanging SmartGrid weighing pan, automatic draft-shield doors, FACT internal adjustment, a 4.5-inch color touchscreen, and LabX compatibility for documented methods and results.

Key specifications

Maximum capacity: 120 g

Readability: 0.1 mg

Settling time: 1.5 s (typical repeatability 0.04 mg; USP minimum weight 82 mg typical)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com