

## Microtrac S3500



Microtrac S3500 Particle Analyzer

Rating: Not Rated Yet

**Price**

Sales price \$22,500.00

[Ask a question about this product](#)

Manufacturer [Microtrac](#)

Description

**!!! SOLD !!!**

This Microtrac S3500 particle analyzer has patented Tri-Laser technology, multi-detector optical system, 0.70 to 1000 micron measurement range, for wet and dry samples, w/SDC, TurboTrac, & Ultra Small Volume Recirculator, current model, lightly used, w/PC & software, over \$60K when purchased new in late 2007. Guaranteed installable by manufacturer (purchase service contract from Microtrac separately)

**Please call us at (650) 871-5707 to check availability of this product!**