

## Millipore Labscale TFF System 29751



500 mL reservoir; Min. recirculating volume: 20 mL

Rating: Not Rated Yet

### Price

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [Millipore](#)

### Description

High-performance, tangential flow filtration (TFF) system ideally suited for process development at volumes from 50 to 1000 mL; process scale-up assisted by built-in controls, and sample preparation. The Labscale TFF System can be operated in either concentration or continuous diafiltration modes.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc., has been based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)

