

Miltenyi AutoMACS



This autoMACS™ Separator is a benchtop instrument for high-speed automated cell sorting. Employing the MACS® Technology, the autoMACS Separator is designed for positive selection as well as depletion of magnetically labeled cells.

Rating: Not Rated Yet

Price

Sales price \$9,750.00

[Ask a question about this product](#)

Manufacturer [Miltenyi](#)

Description

!!! SOLD !!!

The autoMACS™ Separator is a benchtop instrument for high-speed automated cell sorting. Employing the MACS® Technology, the autoMACS Separator is designed for positive selection as well as depletion of magnetically labeled cells. The autoMACS Separator is operated with pre-set separation programs, thus, allowing optimization of cell sorting approaches according to cell abundance and the intensity of marker expression. The separated cells are immediately ready for experiments, cell analysis, or further subset sorting. Automated system cleaning allows convenient maintenance. The autoMACS Separator fits in a standard laminar flow hood or safety cabinet. 100-240V. Comes from a working lab.

Please call us at (650) 871-5707 to check availability of this product!