

## MJ Research PTC-200



MJ Research Gradient PTC-200 with a 96-well Alpha Unit is very clean and comes with our 90-day warranty.

Rating: Not Rated Yet

**Price**

Sales price \$1,850.00

[Ask a question about this product](#)

Manufacturer [MJ Research](#)

### Description

This MJ Research PTC-200 thermal cycler with a 96-well Alpha Unit is very clean and comes with our 90-day warranty. It is Gradient enabled, has a heated lid, temp range of -5 to +105°C, ramps up to 3°C/sec & can store up to 200 programs.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

[Biotech Equipment Sales, Inc.](#)  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)