

Molecular Devices Flexstation



Molecular Devices Flexstation scanning microplate fluorometer with integrated fluid transfer workstation. Comes w/PC loaded w/SoftMax Pro software

Rating: Not Rated Yet

Price

Sales price \$14,250.00

[Ask a question about this product](#)

Manufacturer [Molecular Devices](#)

Description

SOLD!!!

This Molecular Devices Flexstation scanning microplate fluorometer with integrated fluid transfer workstation comes w/PC loaded w/SoftMax Pro software. It's capable of conducting endpoint, kinetic, lambda-scan and well-scan experiments in a multi-well plate format (6, 12, 24, 48 and 96-well) and fluids transfer experiments in 96-well format. Does not accept 384-well plates.