

## Neslab RTE-100



Neslab RTE-100 economical analog refrigerating heating circulating bath, 5L, NEW

Rating: Not Rated Yet

**Price**

Sales price \$950.00

[Ask a question about this product](#)

Manufacturer [Neslab](#)

### Description

**!!!! - RECENTLY SOLD - !!!!**

This Neslab RTE-100 economical analog refrigerating heating circulating bath, has a reservoir capacity of 5L, 350W of cooling at 20C, a temp range of -12 to 100C, can be used for both immersion and external applications, older model but BRAND NEW