

## New Brunswick Innova 2300



This New Brunswick Innova 2300 reciprocal platform shaker has a 36" x 24" platform, digital display of time, speed and alarms, and orbits at 25 to 500 rpm.

Rating: Not Rated Yet

**Price**

Sales price \$2,750.00

[Ask a question about this product](#)

Manufacturer [New Brunswick Scientific](#)

Description

**!!!! - RECENTLY SOLD - !!!!**

This New Brunswick Innova 2300 reciprocal platform shaker has a 36" x 24" platform, digital display of time, speed and alarms, and orbits at 25 to 500 rpm.