

New Brunswick Innova 4230 Refrigerated Incubator Shaker



New Brunswick Scientific Innova 4230 Refrigerated Benchtop Incubator

Rating: Not Rated Yet

Price

Sales price \$3,500.00

[Ask a question about this product](#)

Manufacturer [New Brunswick Scientific](#)

Description

SOLD!!!

New Brunswick Scientific Innova 4230 Refrigerated Benchtop Incubator

The New Brunswick Innova 4230 Refrigerated Incubator Shaker features both heating and refrigeration abilities, a heavy gauge stainless steel frame, and a strong internal shaking platform.

Features:

- double-paned glass window for internal sample viewing
- two adjustable stainless steel shelves and a fluorescent light for quality visibility
- internal AC receptacle for the powering of internal devices
- built-in shaker comes with an audible and visible alarm, control, display, and timer
- internal shaker drive is a triple eccentric counter-balanced drive with 9 sealed-for-life ball bearings and has a motion of ¾ circular standard
- rpm range is between 25 and 400 non stacked, and 25-300 stacked
- onboard timer has a time range of up to 99.9 hours setable in 0.1 hours
- temperature of the 4230 ranges from 20°C below ambient with a minimum of 4°C to maximum of 75°C for 50 Hz units or 80°C for 60 Hz units
- temperature alarms provide indication when temperature deviates > 1.0°C