

New Brunswick Innova 43 Incubator Shaker



30 x 18 in platform

Rating: Not Rated Yet

Price

Sales price \$6,250.00

[Ask a question about this product](#)

Manufacturer [New Brunswick Scientific](#)

Description

The Innova 43 Incubator Shaker is the perfect choice for handling large quantities or large volumes with precision and ease. The cast iron triple-eccentric counterbalanced drive is specifically designed to support high-speed applications and heavy workloads. A Proportional/Integral microprocessor controls allow for uniform motion and accurate regulation of shaking speed. Audible and visual alarms come standard with built-in timer.

Platform Area (ft²): 3.75
Platform Dimensions: 30 x 18 in
RPM - Max: 300, 500
RPM - Min: 25
Orbit Diameter (mm): 2.5 mm

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com