

## New Brunswick Scientific Innova 44R Refrigerating Incubator Shaker



Holds 5L flasks.

Rating: Not Rated Yet

**Price**

Sales price \$9,500.00

[Ask a question about this product](#)

Manufacturer [New Brunswick Scientific](#)

**Description**

Model: Innova 44R

Stackable incubator shaker designed to support high-speed applications and heavy workloads, taking flasks up to 5L. The simple user interface allows for easy programming of speed, temperature and time.

**Specifications:**

Load Bearing Capacity (Metric): 15.5 kg

Speed Range: 25 to 400 rpm

Holds: 5 L Flasks

Load Bearing Capacity (English): 34 lb.

Dimensions (platform) L x W: 30 x 18 in.  
Dimensions (exterior): 33.2 x 53.1 x 26 in.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)