

Nuair Inc. Double Stack In-VitroCell ES NU-5800 Incubator



200L; 7.06 cu ft.

Rating: Not Rated Yet

Price

Sales price \$6,750.00

[Ask a question about this product](#)

Manufacturer [Nuair](#)

Description

The In-VitroCell NU-5800 provides the basic functionality to promote tissue cell growth by controlling Temperature and CO₂ gas levels in the growth chamber. Relative humidity up to 95% is produced through convection by a water pan. The Closed Loop HEPA Filtration System creates clean conditions to minimize potential contamination entering the growth chamber.

Specifications:

Capacity: 200L; 7.06 cu ft.

Relative humidity: 95%

Temperature range: 5C above ambient to 55C

Temperature control: direct heat

Sensor type: Infrared (NDIR)

Interior: Stainless steel

Shelves: 4

Electrical: 115V/60Hz

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com