

Ohaus Benchtop Digital Orbital Shaker



Speed range 50–300 rpm

Rating: Not Rated Yet

Price

Sales price \$400.00

[Ask a question about this product](#)

Manufacturer [Ohaus](#)

Description

The Ohaus Digital Orbital Shaker is a compact, front-controlled benchtop shaker designed for reliable mixing of flasks, tubes, microplates, and culture vessels. The digital interface provides precise control of speed and timed operation, while the ribbed rubber load platform keeps containers stable during agitation. Ideal for microbiology, chemistry, and general lab sample agitation.

Specifications

Motion Type: Orbital shaking

Speed Control: Digital, adjustable (typical range 50–300 rpm — varies by configuration)

Platform Style: Ribbed rubber nonslip platform suitable for a wide range of vessels

Power Input: 120V, 1.8A, 50/60 Hz, 200 W

Controls: Digital speed and time control with independent displays and run/pause control

Timer: Programmable run time or continuous operation

Construction: Heavy-duty steel chassis with front digital keypad

Certifications: TUV listed and FCC Part 15 compliant

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax

sales@biotechequipmentsales.com