

Ohaus Compass CX221 Portable Scale



Maximum capacity: 220 g

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Ohaus](#)

Description

The Ohaus Compass CX221 is a compact portable scale for basic weighing in classrooms, labs, and light industrial use. It has a low-profile stackable ABS housing, a removable stainless-steel pan, a backlit LCD, two-button operation, push-button external calibration, and overload protection.

Key specifications

Maximum capacity: 220 g

Readability: 0.1 g

Stabilization time: 1.5 s (repeatability 0.1 g; linearity ± 0.1 g)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com