

OHAUS CS-200 Compact Digital Balance



Capacity: 200 g

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Ohaus](#)

Description

The OHAUS CS-200 is a lightweight, portable digital balance designed for general-purpose weighing in laboratory, educational, and industrial settings. It features a stainless steel weighing pan, LCD display, and battery or AC power operation.

Key specifications

Capacity: 200 g

Readability: 0.1 g

Pan dimensions: 145 × 133 mm (5.75 × 5.25 in.)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com