

Ohaus Starter 2200 Bench Meter - NEW in Box



User-friendly.

Rating: Not Rated Yet

Price

Sales price \$200.00

[Ask a question about this product](#)

Manufacturer [Ohaus](#)

Description

Designed with a clear, user-friendly interface and the ability to store up to 99 items in its internal memory, the Starter 2200 pH bench meter is a reliable and efficient addition to basic laboratories or science classrooms. An affordable, intuitive pH measurement instrument, the Starter 2200 pH bench meter also features auto-buffer recognition & auto-endpoint modes that make measurement and calibration simple, straightforward and efficient.

Specifications:

Model: ST2200

pH
pH Measuring Range: 0.00 to 14.00
pH Resolution: 0.01 pH
Selectable Resolution: No
Accuracy: ± 0.01 pH
Pre-Defined Buffer Groups: 1

ORP, Redox
Measuring Range: ± 1999 mV
Resolution 1 mV

Meters: Ohaus Starter 2200 Bench Meter - NEW in Box

Accuracy: ± 1.5 mV
Units mV

Temperature
Measuring Range: 0.0 to 100.0 °C, 32.0 °F to 212.0 °F
Resolution: 0.1 °C, 0.1 °F
Accuracy: ± 0.5 °C, ± 0.5 °F

Calibration
Calibration: NoCalibration
Calibration points: 1 to 3 points
Calibration: Sign Face Sign
Calibration mode: Linear

Meter Specifications
Display Type: 3.5" segment LCD with backlight
Measurement End-point Modes : Auto-stop, Continuous
Datalog for measurement: 99 sets
Datalog for calibration: Last calibration
Keypad: Membrane
pH electrode input: BNC
Temperature input: Cinch, NTC 30 k Ω
Installation Overvoltage: Category II
Pollution: Degree 2
Power Supply: Universal, 100-240 VAC, 50-60 Hz
Ambient Operating: 5 to 40 °C, 5 to 80 %, non-condensing

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com