


Omni International Beadruptor 12



Process up to 12 samples at one time.

Rating: Not Rated Yet

Price

Sales price \$3,150.00

[Ask a question about this product](#)

Description

- Process up to 12 samples at one time
- Non-stop processing with no cool-down between runs
- No cross-contamination
- Small footprint

Powerful, fast, and economical, the Omni Bead Ruptor 12 is ideal for many applications that require grinding, lysing, or homogenization of biological samples prior to molecular extraction. Convenient disposable sample tubes come pre-filled with a variety of lysing beads. Up to twelve 2 mL and/or four 7 mL tubes can be processed at once for efficient and consistent high-yield homogenization. Unlike many bead mill homogenizers, the Bead Ruptor 12 does not require a cool-down period between runs. This means more samples can be processed in a smaller amount of time. Tubes are easily loaded into the angled, front-facing holder for optimal bead/sample interactions.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com