

Panasonic MLS-3781L Autoclave



230V; 75L working capacity

Rating: Not Rated Yet

Price

Sales price \$9,250.00

[Ask a question about this product](#)

Manufacturer [Panasonic](#)

Description

Model: MLS-3781L-PA

Specifications:

Power supply: 230V, 50/60 Hz, 17.4 A

Power consumption: 4 kW

Autoclaves: Panasonic MLS-3781L Autoclave

External dimensions: 47B(W) x 632(D) x 965(H) mm
Internal capacity/ effective: 80L; 75L
Chamber dimensions: 370 (dia.) x630(D); (Effective depth including cover: 6BB) mm
Chamber material: Austenitic stainless steel
Sterilization temperature: 115°C - 135°C
Melting temperature: 60°C - 114°C
Heat-retention Temperature: 45°C - 60°C
Thermometer: Thermistor, digital display (25°C - 141°C)
Safety Valve releasing pressure: 240kPa (34.Bpsi, 2.4bar)
Pressure gauge range: -0.1 - 0.4 MPa
Timer: Sterilization and Melting - 1 min ~ 5 hours; Heat retent: 72 hours, fixed; Timer: 1 week

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com