


## Perkin Elmer OmniPrep 96 Automated Homogenizer with Tablet



Variable speed from 500 – 28,000 rpm

Rating: Not Rated Yet

**Price**

Sales price \$63,995.00

[Ask a question about this product](#)

Manufacturer [Perkin Elmer](#)

### Description

This fully automated homogenization workstation enables walk-away processing in your lab. The small footprint fits into most biosafety cabinets and fume hoods. User-friendly programming makes saving, recalling, or changing existing methods easy. Omni patented Whisper Drive Technology with brushless motors delivers maintenance-free, ultra-quiet, maximum processing power.

Part number: 51-02A-1

### Features:

- Automated homogenization of up to ninety-six samples
- Process 250  $\mu$ L to 40 mL in 2 to 50 mL tubes
- Create custom sample processing methods in the Prep 96 software
- Unique oscillating motion to ensure complete homogenization
- Variable speed from 500 – 28,000 rpm

- Compatible with 7 and 12 mm OMNI Tips Plastic Disposable Probes
- Compatible with 5, 7, 10 mm Stainless Steel Generator Probes
- Compatible with 7 and 12 mm OMNI Tip Hybrid Probes

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)