

## PHCBI MIR-554-PA: MIR Series 14.3 cu.ft. Large Capacity Refrigerating Incubator -- New!



Brand new! Capacity: 14.3 cu.ft. (406 L)

Rating: Not Rated Yet

**Price**

Sales price \$5,950.00

[Ask a question about this product](#)

### Description

Uniform and consistent temperature control through an easy-to-use, programmable interface for a wide range of heated and cooled applications. Stepwise, temperature and lighting programming capabilities for extensive protocols in a wide range of applications. Forced air circulation.

### Specifications:

Capacity: 14.3 cu.ft. (406 L)

Exterior Dimensions: 31.5" W x 32.8" D x 71.3" H (800 x 832 x 1810 mm)

Interior Dimensions: 25.2" W x 21.7" D x 45.7" H (640 x 550 x 1160 mm)

Shelves: As pictured

Temperature control range: -10°C to 60°C (ambient temperature: 5°C to 35°C, no load)

Temperature uniformity:  $\pm 0.5^{\circ}\text{C}$  (SV:37%, ambient temperature: 20°C, no load)

Alarms: High/Low temp alarm, independent over-heat/over-cool limiter

Net weight: 430 lbs (195 kg)

Electrical requirements: 290 W, AC 115 V, 60 Hz

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-

quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)